



# DuCanto Catheter™

## SSCOR DuCanto Catheter/Yankauer Comparison



### Standard Yankauer

- Designed in 1907 for removing blood and saliva during tonsillectomy
- Small distal openings to permit fluid removal
- Often has a thumb port so that suction can be applied and interrupted

### SSCOR DuCanto Catheter

- Designed in 2016 to remove material from the upper airway during resuscitation
- Large lumen (.26") for removal of emesis, clots, tissue during resuscitation
- No thumb port so that suction is 'always on' and removing material from the airway even when the caregiver's hands are occupied with other aspects of the resuscitation
- Anatomical shape allows the distal tip of the SSCOR DuCanto Catheter to rest near the top of the esophagus to continually drain the hypopharynx during intubation

 DuCanto Catheter™



Clearing The Airway Is Our #1 Priority

[sscor.com](http://sscor.com)

SSCOR, Inc.  
2125 Madera Rd Simi Valley, CA 93065  
email: [sales@sscor.com](mailto:sales@sscor.com) or [techsupport@sscor.com](mailto:techsupport@sscor.com)  
+1 818 504 4054  
[www.sscor.com](http://www.sscor.com)



LIT-200-00002-1 Rev B